Data sheet

Butyl SepFast[™] (4HF, 6HF, Large Beads)

Phenyl SepFast[™] (HighRes, 4HF, 6HF, Large Beads)

Phenyl SepFast[™] 6HF High sub

Octyl SepFast[™] (4HF, 6HF, Large Beads)

1. Introduction

Butyl, Phenyl and Octyl SepFast are hydrophobic interaction chromatography (HIC) adsorbents having a carbon chain of C4, benzene ring and C8, respectively. They are specially designed for the purification of biological molecules based on their hydrophobicity profiles.

HIC is a versatile technique and could show high selectivity to individual molecules according to their exposed hydrophobic zones. It is particularly useful for intermediate and final-stage purifications. A HIC medium normally binds at moderate to high salt concentrations. It is logical to place HIC step after an IEX step where molecules are usually eluted at high salt conditions.

HIC media shows much milder purification conditions than reversed phase chromatography (RPC) media. Better biological activity could be maintained in HIC operations than RPC operations.

SepFast HIC media offers a broad range of choices. Generally speaking, the longer the carbon chain, the higher the surface hydrophobicity. The base matrix is made of agarose that has been highly cross-linked. It is very stable to most of the chemical conditions experienced in bioprocessing industry. Process developers and manufacturers have the chance to choose the best HIC media for their applications. The core advantages are:

- High sample loading capacity
- More choices than other suppliers
- High separation power
- All versions can be scaled up and produced for GMP use
- Regulatory support files available

The feature and selection guide are listed below:

4HF serial	6HF serial	Large Beads	HighRes
(50 – 150 μm)	(50 – 150 μm)	(150 – 350 μm)	(25 – 50 μm)
The above HIC media is designed to purify most medium to large proteins	The above HIC media is designed to purify peptides or smaller proteins	The above HIC media is designed to purify proteins from crude or viscous samples	The above HIC media is designed to purify proteins that require high resolution

Characteristics of SepFast HIC media:

	Butyl	Phenyl	Octyl
Matrix	Highly cross-linked agarose; 4% agarose for SepFast 4HF, 6% agarose for SepFast 6HF		
Functional group	-(CH ₂) ₃ CH ₃	-C ₆ H ₆	-(CH ₂) ₇ CH ₃
Ligand density	Around 40 μmol/ml	Around 25 μmol/ml Around 40 μmol/ml (high sub)	Around 4 μmol/ml
Binding capacity	>20 mg lyzosyme / ml	>20 mg BSA / ml	>20 mg BSA / ml
Particle size	25 - 50 μm (HighRes), 50 - 150 μm (4HF and 6HF), 150 - 350 μm (Large Beads)		
Pressure-flow property*	>100 cm/h for HighRes; >300 cm/h for 4HF serials; >500 cm/h for 6HF serials; >1500 cm/h for Large Beads		
Operational pressure	Up to 3 bar		
pH stability	2-14 (short term) and 3-12 (long term)		
Working temperature	+4°C to +30°C		
Chemical stability	All commonly used buffers; 1 M acetic acid, 1 M NaOH, 6M guanidine hydrochloride, 8 M urea, 30% acetonitrile, 30% isopropanol, 70% ethanol, 3 M (NH ₄) ₂ SO ₄		
Storage	20% ethanol		

^{*}Measured in a 32 mm ID column at a bed height of 15 cm.

2. Method optimization

We recommend scouting the parameters among loading capacity, flow velocity, binding pH, binding ionic strength, elution speed and gradient etc. We recommend paying special attention to optimize elution conditions to achieve the best separation power.

In general, balancing the degree of component separation against process throughput is the major consideration when optimizing a method. Besides, for the purification of instable or shearing-force sensitive molecules, the operational condition needs be optimised to balance the throughput and the possible damage to the target molecule.

3. Maintenance

Depending on the individual applications, the media may be used many times. For the re-use purpose, please see the following instructions.

Cleaning-in-place (CIP)

CIP is a procedure that removes strongly bound materials such as lipids, endotoxins and denatured proteins that remain in the adsorbent surface after the regeneration. Regular CIP prevents the build up of contaminants in the packed bed and helps to maintain the column performance.

A specific CIP protocol should be developed for each process according to the type of contaminants present. The frequency of CIP depends on the nature of individual applications.

The following information works as a general guidance.

The contaminants bound by hydrophobic nature can be removed by the following reagents: 1 M NaOH, low percentage non-ionic detergents (e.g. 0.1 - 2%), 30% isopropanol in basic or acidic conditions (e.g. in the presence of acetic acid or phosphoric acid). A combination of the above reagents can be explored as well. In general, the incubation time should be longer (e.g. from 30 minutes to 2 hours) to ensure full dissociation of the contaminants.

Sanitization

Sanitization using 0.5-1.0 M NaOH with a contact time of 30 mins is recommended.

4. Storage

The media should be stored in 20% ethanol or 0.02% sodium azide to prevent microbial growth. Store the media at a temperature of +4°C to +30°C. Before use, equilibrate the media with at least 5 bed volumes of the running buffer.

5. Ordering information

Butyl 6HF, Phenyl 4HF, Phenyl 6HF, Phenyl 6HF, Phenyl 6HF, Phenyl 6HF high sub, Octyl 4HF and Octyl 6HF, respectively	Product	Quantity	Code no.
100 ml	HIC SepFast Selection Kit	Butyl 6HF, Phenyl 4HF, Phenyl 6HF, Phenyl 6HF high sub, Octyl 4HF and Octyl 6HF,	480000-7x1ML
1 litre	Butyl SepFast 4HF	25 ml	480101-25ML
Pre-packed column 5 x 1 ml 5 x 5 ml 480101-5x1ML 480101-5x5ML Butyl SepFast 6HF 25 ml 480301-25ML 480301-25ML 480301-10ML 1 litre Pre-packed column 5 x 1 ml 480301-5x1ML 480301-5x1ML 480301-5x1ML 5 x 5 ml Butyl SepFast Large Beads 25 ml 480201-5x5ML 480201-25ML 480201-10ML 1 litre Phenyl SepFast 4HF 25 ml 481001-25ML 480201-10ML 1 litre Pre-packed column 5 x 1 ml 481001-10ML 1 litre Pre-packed column 5 x 1 ml 481001-5x1ML 1 litre Phenyl SepFast 6HF 25 ml 481201-5x1ML 1 litre Phenyl SepFast 6HF 25 ml 481201-25ML 1 litre Pre-packed column 5 x 1 ml 481201-10ML 1 litre Pre-packed column 5 x 1 ml 481201-5x1ML 1 litre Pre-packed column 5 x 1 ml 481201-5x1ML 1 litre Pre-packed column 5 x 1 ml 481201-5x1ML 1 litre Pre-packed column 5 x 1 ml 481201-5x1ML 1 litre Pre-packed column 5 x 1 ml 481201-5x1ML 1 litre Pre-packed column 5 x 1 ml 481201-5x1ML 1 litre Pre-packed column 5 x 1 ml		100 ml	480101-100ML
S x 5 ml		1 litre	480101-1L
Butyl SepFast 6HF 25 ml	Pre-packed column	5 x 1 ml	480101-5x1ML
100 ml		5 x 5 ml	480101-5x5ML
1 litre	Butyl SepFast 6HF	25 ml	480301-25ML
Pre-packed column 5 x 1 ml 480301-5x1ML 5 x 5 ml 480301-5x5ML Butyl SepFast Large Beads 25 ml 480201-25ML 100 ml 480201-100ML 1 litre 480201-1L Phenyl SepFast 4HF 25 ml 481001-25ML Pre-packed column 1 litre 481001-1L Pre-packed column 5 x 1 ml 481001-5x1ML 5 x 5 ml 481001-5x5ML Phenyl SepFast 6HF 25 ml 481201-25ML Pre-packed column 5 x 1 ml 481201-1L Pre-packed column 5 x 1 ml 481201-1L Pre-packed column 5 x 1 ml 481201-5x1ML Phenyl SepFast 6HF High sub 25 ml 481202-25ML Phenyl SepFast 6HF High sub 25 ml 481202-25ML 100 ml 481202-100ML 481202-100ML		100 ml	480301-100ML
S x 5 ml		1 litre	480301-1L
Butyl SepFast Large Beads 25 ml	Pre-packed column	5 x 1 ml	480301-5x1ML
100 ml		5 x 5 ml	480301-5x5ML
1 litre	Butyl SepFast Large Beads	25 ml	480201-25ML
Phenyl SepFast 4HF 25 ml 481001-25ML 100 ml 481001-100ML 1 litre 481001-5x1ML 5 x 1 ml 5 x 5 ml 481001-5x5ML Phenyl SepFast 6HF 25 ml 481201-25ML 100 ml 481201-100ML 1 litre 481201-11L Pre-packed column 5 x 1 ml 481201-11L Pre-packed column 5 x 1 ml 481201-5x5ML Phenyl SepFast 6HF High sub 25 ml 481201-5x5ML Phenyl SepFast 6HF High sub 25 ml 481202-25ML 481202-100ML		100 ml	480201-100ML
100 ml 481001-100ML 1 litre 481001-1L Pre-packed column 5 x 1 ml 481001-5x1ML 5 x 5 ml 481001-5x5ML Phenyl SepFast 6HF 25 ml 481201-25ML 100 ml 481201-100ML 1 litre 481201-1L Pre-packed column 5 x 1 ml 481201-5x1ML 5 x 5 ml 481201-5x1ML 5 x 5 ml 481201-5x5ML Phenyl SepFast 6HF High sub 25 ml 481202-25ML 100 ml 481202-25ML		1 litre	480201-1L
1 litre 481001-1L Pre-packed column 5 x 1 ml 481001-5x1ML 5 x 5 ml 481001-5x5ML Phenyl SepFast 6HF 25 ml 481201-25ML 100 ml 481201-100ML 1 litre 481201-1L Pre-packed column 5 x 1 ml 481201-5x1ML 5 x 5 ml 481201-5x5ML Phenyl SepFast 6HF High sub 25 ml 481202-25ML 100 ml 481202-100ML	Phenyl SepFast 4HF	25 ml	481001-25ML
Pre-packed column 5 x 1 ml 481001-5x1ML 5 x 5 ml 481001-5x5ML Phenyl SepFast 6HF 25 ml 481201-25ML 100 ml 481201-100ML 1 litre 481201-1L Pre-packed column 5 x 1 ml 481201-5x1ML 5 x 5 ml 481201-5x5ML Phenyl SepFast 6HF High sub 25 ml 481202-25ML 100 ml 481202-100ML		100 ml	481001-100ML
5 x 5 ml		1 litre	481001-1L
Phenyl SepFast 6HF 25 ml 481201-25ML 100 ml 481201-100ML 1 litre 481201-1L Pre-packed column 5 x 1 ml 481201-5x1ML 5 x 5 ml 481201-5x5ML 481201-5x5ML Phenyl SepFast 6HF High sub 25 ml 481202-25ML 100 ml 481202-100ML	Pre-packed column	5 x 1 ml	481001-5x1ML
100 ml 481201-100ML 1 litre 481201-1L Pre-packed column 5 x 1 ml 481201-5x1ML 5 x 5 ml 481201-5x5ML Phenyl SepFast 6HF High sub 25 ml 481202-25ML 100 ml 481202-100ML		5 x 5 ml	481001-5x5ML
$\begin{array}{cccccccccccccccccccccccccccccccccccc$	Phenyl SepFast 6HF	25 ml	481201-25ML
Pre-packed column 5 x 1 ml 481201-5x1ML 5 x 5 ml 481201-5x5ML Phenyl SepFast 6HF High sub 25 ml 481202-25ML 100 ml 481202-100ML		100 ml	481201-100ML
5 x 5 ml 481201-5x5ML Phenyl SepFast 6HF High sub 25 ml 481202-25ML 100 ml 481202-100ML		1 litre	481201-1L
Phenyl SepFast 6HF High sub 25 ml 481202-25ML 100 ml 481202-100ML	Pre-packed column	5 x 1 ml	481201-5x1ML
100 ml 481202-100ML		5 x 5 ml	481201-5x5ML
	Phenyl SepFast 6HF High sub	25 ml	481202-25ML
1 litre 481202-1L		100 ml	481202-100ML
		1 litre	481202-1L

Pre-packed column	5 x 1 ml	481202-5x1ML
	5 x 5 ml	481202-5x5ML
Phenyl SepFast HighRes	25 ml	481101-25ML
	100 ml	481101-100ML
	1 litre	481101-1L
Pre-packed column	5 x 1 ml	481101-5x1ML
	5 x 5 ml	481101-5x5ML
Phenyl SepFast Large Beads	25 ml	481301-25ML
	100 ml	481301-100ML
	1 litre	481301-1L
Octyl SepFast 4HF	25 ml	481601-25ML
	100 ml	481601-100ML
	1 litre	481601-1L
Pre-packed column	5 x 1 ml	481601-5x1ML
	5 x 5 ml	481601-5x5ML
Octyl SepFast 6HF	25 ml	481801-25ML
	100 ml	481801-100ML
	1 litre	481801-1L
Pre-packed column	5 x 1 ml	481801-5x1ML
	5 x 5 ml	481801-5x5ML
Octyl SepFast Large Beads	25 ml	481701-25ML
	100 ml	481701-100ML
	1 litre	481701-1L



BioToolomics Ltd Unit 30A, Number 1 Industrial Estate Consett County Durham, DH8 6TJ United Kingdom

www.biotoolomics.com

Registered or registration-pending trademark of BioToolomics Litd: BioToolomics, SepFast.

All goods and services are sold subject to the terms and conditions of sale of BioToolomics Ltd. The company reserves the rights, subject to regulatory or contractual approval, if required, to make changes in the specifications and features shown herein, or discontinue the products described at any time without notice or obligation. Contact BioToolomics Ltd for the most current information.