

Stop Plug



Stop Plug has 10-32 UNF male thread. For storage purpose, it can be used for the fingertight seal of packed columns. The Stop Plug can be used in our 6.2 mm i.d. column, 11 mm i.d. column and 16 mm i.d. column. It is compatible to GE Healthcare small columns.

BioToolomics Ltd
Unit 30A, Number 1 Industrial Estate
Consett
County Durham, DH8 6TJ
United Kingdom

Registered or registration-pending trademark of BioToolomics Ltd: BioToolomics, SepFast.

All goods and services are sold subject to the terms and conditions of sale of BioToolomics Ltd. The company reserves the rights, subject to regulatory or contractual approval, if required, to make changes in the specifications and features shown herein, or discontinue the products described at any time without notice or obligation. Contact BioToolomics Ltd for the most current information.

© 2008-2012 BioToolomics Ltd – All rights reserved.